

ESM Table 2. Sequences of custom-designed primers (forward and reverse) and probes used for quantitative real-time PCR

Gene		Sequence (5'-3')
<i>Cycs</i>	Forward	GGACCAAATCTCCACGGTCT
	Reverse	CCCCAGGTGATGCCTTTGT
	Probe	AGCCTGGCCTGTCTTCCGCCC
<i>Fasn</i>	Forward	CTTAGCAGAGATCCCGAGACG
	Reverse	GGTCCTTTGAAGTCGAAGAAGAA
	Probe	TGGGCTACAGCATGGTGGGCTG